




FGR Antibody

Product Code	CSB-PA824952
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09769
Immunogen	Synthesized peptide derived from internal of Human FGR.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total FGR protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	C-FGR; EC 2.7.10.2; FGR protein; P55-FGR; Proto-oncogene tyrosine-protein kinase FGR
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FGR

Image



Western blot analysis of extracts from HUVEC cells and COLO205 cells, using FGR antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.