




LYL1 Antibody

Product Code	CSB-PA284084
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12980
Immunogen	Synthesized peptide derived from internal of Human Lyl-1.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total Lyl-1 protein.
Tested Applications	ELISA, IHC; IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Protein lyl-1; Lymphoblastic leukemia-derived sequence 1; LYL1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LYL1

Image



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using Lyl-1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.