





TPT1 Antibody

Product Code	CSB-PA588140
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13693
Immunogen	Synthesized peptide derived from internal of Human TCTP.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total TCTP protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	HRF; TCTP; fortilin; histamine-releasing factor; translationally-controlled tumor protein
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TPT1

Image



Western blot analysis of extracts from HepG2 cells, COLO205 cells and HeLa cells, using TCTP antibody.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using TCTP antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.