





IRF2 Antibody

Product Code	CSB-PA072893
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14316
Immunogen	Synthesized peptide derived from internal of Human IRF2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total IRF2 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	IRF2; Interferon regulatory factor 2;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IRF2

Image



Western blot analysis of extracts from Jurkat cells and COLO205 cells, using IRF2 antibody.



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue using IRF2 antibody.