




# HOXB1/HOXA1/HOXD1 Antibody

<b>Product Code</b>	CSB-PA028363
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P14653
<b>Immunogen</b>	Synthesized peptide derived from internal of Human HOXA1/B1/D1.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse
<b>Specificity</b>	The antibody detects endogenous levels of total HOXA1/B1/D1 protein.
<b>Tested Applications</b>	ELISA,IHC;IHC:1:50-1:100
<b>Form</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Homeobox protein Hox-B1; Hox-2l; HOXB1; HOX2l;
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	HOXB1/HOXA1/HOXD1

## Image



Immunohistochemistry analysis of paraffin-embedded human testis tissue, using HOXA1/B1/D1 antibody.