




HOXB1/HOXA1/HOXD1 Antibody

Product Code	CSB-PA028363
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14653
Immunogen	Synthesized peptide derived from internal of Human HOXA1/B1/D1.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total HOXA1/B1/D1 protein.
Tested Applications	ELISA, IHC; IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Homeobox protein Hox-B1; Hox-2l; HOXB1; HOX2l;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HOXB1/HOXA1/HOXD1

Image



Immunohistochemistry analysis of paraffin-embedded human testis tissue, using HOXA1/B1/D1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.