





TCF3 Antibody

Product Code	CSB-PA079787
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15923
Immunogen	Synthesized peptide derived from C-terminal of Human TCF3.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total TCF3 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	E12; ITF1; Immunoglobulin enhancer binding factor E12/E47; Immunoglobulin transcription factor-1; Kappa-E2-binding factor
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TCF3

Image



Western blot analysis of extracts from 293 cells, using TCF3 antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using TCF3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.