




EPCAM Antibody

Product Code	CSB-PA787369
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16422
Immunogen	Synthesized peptide derived from internal of Human TACD1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TACD1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Tumor-associated calcium signal transducer 1 [Precursor]; Major gastrointestinal tumor-associated protein GA733-2; Epithelial cell surface antigen; Epithelial glycoprotein; EGP
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EPCAM

Image



Western blot analysis of extracts from Jurkat cells, HuvEc cells and COLO205 cells, using TACD1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.