




ZNF18 Antibody

Product Code	CSB-PA598851
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17022
Immunogen	Synthesized peptide derived from internal of Human ZNF18.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Specificity	The antibody detects endogenous levels of total ZNF18 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Zinc finger protein 18; Zinc finger protein KOX11; Heart development-specific gene 1 protein; Zinc finger protein with KRAB and SCAN domains 6; ZNF18
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZNF18

Image



Western blot analysis of extracts from HeLa cells, using ZNF18 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.