





UBTF Antibody

Product Code	CSB-PA829201
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17480
Immunogen	Synthesized peptide derived from internal of Human UBF1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total UBF1 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Autoantigen NOR-90; Nucleolar transcription factor 1; TCFUBF; UBF-1; UBF1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	UBTF

Image



Western blot analysis of extracts from HeLa cells and Jurkat cells, using UBF1 antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using UBF1 antibody.