




ADRA2B Antibody

Product Code	CSB-PA585658
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18089
Immunogen	Synthesized peptide derived from internal of Human Adrenergic Receptor alpha-2B.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total Adrenergic Receptor alpha-2B protein.
Tested Applications	ELISA,IF;IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ADA2B; ADRA2L1; ADRA2RL1; ADRARL1; ALPHA2BAR
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ADRA2B

Image



Immunofluorescence analysis of HepG2 cells, using Adrenergic Receptor alpha-2B antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.