




EGR1/EGR2 Antibody

Product Code	CSB-PA253458
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18146
Immunogen	Synthesized peptide derived from internal of Human EGR1/2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total EGR1/2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	AT225; EGR-1; Early growth response protein 1; KROX-24; Krox-24 protein
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EGR1/EGR2

Image



Western blot analysis of extracts from HeLa/Jurkat cells, using EGR1/2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.