





TFE3 Antibody

Product Code	CSB-PA827304
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19532
Immunogen	Synthesized peptide derived from internal of Human TFE3.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total TFE3 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Transcription factor E3; TFE3;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TFE3

Image



Western blot analysis of extracts from K562 cells and A549 cells, using TFE3 antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using TFE3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.