




TBP Antibody

Product Code	CSB-PA276652
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20226
Immunogen	Synthesized peptide derived from internal of Human TBP.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total TBP protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	TATA binding factor; TATA sequence-binding protein; TATA-box binding protein; TATA-box factor; TFIID
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TBP

Image



Western blot analysis of extracts from 293 cells, using TBP antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.