





TNFAIP3 Antibody

Product Code	CSB-PA257212
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21580
Immunogen	Synthesized peptide derived from internal of Human TNAP3.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total TNAP3 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	A20; DNA binding protein A20; EC 3.-.-.-; TNAP3; TNFA1P2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TNFAIP3

Image



Western blot analysis of extracts from RAW264.7 cells, using TNAP3 antibody.



Immunofluorescence analysis of HepG2 cells, using TNAP3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.