




PRKACG Antibody

Product Code	CSB-PA925968
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22612
Immunogen	Synthesized peptide derived from N-terminal of Human KAPCG.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total KAPCG protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	cAMP-dependent protein kinase; gamma-catalytic subunit; EC 2.7.11.11; PKA C-gamma; PRKACG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PRKACG

Image



Western blot analysis of extracts from COLO205 cells, using KAPCG antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.