




PRKACB Antibody

Product Code	CSB-PA987088
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22694
Immunogen	Synthesized peptide derived from internal of Human KAPCB.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total KAPCB protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.11.11; KAPB; KAPCB; PKA C-beta; PKA-beta
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PRKACB

Image



Western blot analysis of extracts from Jurkat cells and COLO205 cells, using KAPCB antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.