





JAK1 Antibody

Product Code	CSB-PA279447
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23458
Immunogen	Synthesized peptide derived from internal of Human JAK1.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total JAK1 protein.
Tested Applications	ELISA, IHC, IF; IHC: 1:50-1:100, IF: 1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.10.2; JAK-1; JAK1; JAK1A; Janus kinase 1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	JAK1

Image



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using JAK1 antibody.



Immunofluorescence analysis of HeLa cells, using JAK1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.