




NFYA Antibody

Product Code	CSB-PA288741
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23511
Immunogen	Synthesized peptide derived from internal of Human NFYA.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total NFYA protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CAAT-box DNA binding protein subunit A; CBF-A; CBF-B; CCAAT-binding transcription factor subunit B; HAP2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NFYA

Image



Western blot analysis of extracts from 293 cells, HepG2 cells and COLO205 cells, using NFYA antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.