





PTX3 Antibody

Product Code	CSB-PA962105
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26022
Immunogen	Synthesized peptide derived from N-terminal of Human PTX3.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total PTX3 protein.
Tested Applications	ELISA, WB, IF; WB: 1:500-1:3000, IF: 1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Pentraxin-related protein PTX3 [Precursor]; Pentaxin-related protein PTX3; Tumor necrosis factor-inducible protein TSG-14; PTX3; TNFAIP5
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PTX3

Image



Western blot analysis of extracts from 293 cells, using PTX3 antibody.



Immunofluorescence analysis of HUVEC cells, using PTX3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.