




MAPK3 Antibody

Product Code	CSB-PA272447
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27361
Immunogen	Synthesized peptide derived from internal of Human p44 MAPK.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total p44 MAPK protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.11.24; ERK-1; ERK1; ERT2; Extracellular signal- regulated kinase 1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MAPK3

Image



Western blot analysis of extracts from Jurkat cells and HUVEC cells, using p44 MAPK antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.