





MARK3 Antibody

Product Code	CSB-PA052531
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27448
Immunogen	Synthesized peptide derived from N-terminal of Human MARK3.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total MARK3 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	C-TAK1; Cdc25C- associated protein kinase 1; EC 2.7.11.1; MAP/microtubule affinity-regulating kinase 3; Protein kinase STK10
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MARK3

Image



Western blot analysis of extracts from HeLa cells and 293 cells, using MARK3 antibody.



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using MARK3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.