




EPHA3 Antibody

Product Code	CSB-PA298527
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29320
Immunogen	Synthesized peptide derived from internal of Human EPHA3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total EPHA3 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EPA3; ETK; ETK1; Ephrin type-A receptor 3 precursor; HEK
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EPHA3

Image



Western blot analysis of extracts from 293 cells, using EPHA3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.