




LTK Antibody

Product Code	CSB-PA285877
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29376
Immunogen	Synthesized peptide derived from internal of Human LTK.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total LTK protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.1.112; EC 2.7.10.1; Leukocyte tyrosine kinase receptor precursor; Protein tyrosine kinase-1; TYK1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LTK

Image



Western blot analysis of extracts from Jurkat cells, using LTK antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.