





PPIF Antibody

Product Code	CSB-PA929434
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30405
Immunogen	Synthesized peptide derived from internal of Human PPIF.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total PPIF protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Peptidyl-prolyl cis-trans isomerase; mitochondrial [Precursor]; EC 5.2.1.8; PPIase; Rotamase
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PPIF

Image



Western blot analysis of extracts from HeLa cells, 293 cells and HepG2 cells, using PPIF antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using PPIF antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.