





PPIF Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA929434 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P30405 |
| Immunogen | Synthesized peptide derived from internal of Human PPIF. |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Specificity | The antibody detects endogenous levels of total PPIF protein. |
| Tested Applications | ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100 |
| Form | Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Alias | Peptidyl-prolyl cis-trans isomerase; mitochondrial [Precursor]; EC 5.2.1.8; PPIase; Rotamase |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | PPIF |

Image



Western blot analysis of extracts from HeLa cells, 293 cells and HepG2 cells, using PPIF antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using PPIF antibody.