




PRKAR1B Antibody

Product Code	CSB-PA129156
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31321
Immunogen	Synthesized peptide derived from internal of Human KAP1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total KAP1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	KAP1; PRKAR1B; cAMP-dependent protein kinase type I-beta regulatory subunit;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PRKAR1B

Image



Western blot analysis of extracts from COLO205 cells, using KAP1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.