





GTF2H1 Antibody

Product Code	CSB-PA934028
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32780
Immunogen	Synthesized peptide derived from N-terminal of Human TF2H1.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total TF2H1 protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:3000, IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	TFIIH basal transcription factor complex p62 subunit; Basic transcription factor 62 kDa subunit; BTF2-p62; General transcription factor IIH polypeptide 1; GTF2H1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GTF2H1

Image



Western blot analysis of extracts from Jurkat cells, using TF2H1 antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using TF2H1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.