




TNFSF8 Antibody

Product Code	CSB-PA951710
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32971
Immunogen	Synthesized peptide derived from internal of Human CD153.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total CD153 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Tumor necrosis factor ligand superfamily member 8; CD30 ligand; CD30-L; CD153 antigen; TNFSF8
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TNFSF8

Image



Western blot analysis of extracts from RAW264.7 cells, using CD153 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.