




GRK5 Antibody

Product Code	CSB-PA794868
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P34947
Immunogen	Synthesized peptide derived from internal of Human GRK5.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total GRK5 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.11.16; G-protein-coupled receptor kinase 5; GPRK5; kinase GRK5;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GRK5

Image



Western blot analysis of extracts from HUVEC cells, HepG2 cells and HeLa cells, using GRK5 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.