




HOXD12 Antibody

Product Code	CSB-PA828729
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35452
Immunogen	Synthesized peptide derived from internal of Human HOXD12.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total HOXD12 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Homeobox protein Hox-D12; Hox-4H; HOXD12; HOX4H;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HOXD12

Image



Western blot analysis of extracts from K562, JK and 3T3 cells using HOXD12 antibody.