





CSK Antibody

Product Code	CSB-PA283003
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41240
Immunogen	Synthesized peptide derived from internal of Human CSK.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total CSK protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Tyrosine-protein kinase CSK; EC 2.7.10.2; C-SRC kinase; Protein-tyrosine kinase CYL; CSK
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CSK

Image



Western blot analysis of extracts from Jurkat cells, using CSK antibody.



Immunofluorescence analysis of A549 cells, using CSK antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.