





TNFSF9 Antibody

Product Code	CSB-PA203832
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41273
Immunogen	Synthesized peptide derived from internal of Human TNFSF9.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TNFSF9 protein.
Tested Applications	ELISA, WB, IF; WB: 1:500-1:3000, IF: 1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Tumor necrosis factor ligand superfamily member 9; 4-1BB ligand; 4-1BBL; TNFSF9;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TNFSF9

Image



Western blot analysis of extracts from HUVEC cells, using TNFSF9 antibody.



Immunofluorescence analysis of HeLa cells,
using TNFSF9 antibody.