




MAP2K4 Antibody

Product Code	CSB-PA291313
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45985
Immunogen	Synthesized peptide derived from internal of Human JNKK.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total JNKK protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	C-JUN N-terminal kinase kinase 1; Dual specificity mitogen-activated protein kinase kinase 4; EC 2.7.12.2; JNK activating kinase 1; JNK kinase 1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MAP2K4

Image



Western blot analysis of extracts from Jurkat cells and 293 cells, using JNKK antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.