




PHKA1 Antibody

Product Code	CSB-PA566650
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46020
Immunogen	Synthesized peptide derived from internal of Human KPBI/2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total KPBI/2 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	KPB1; PHKA; Phosphorylase B kinase alpha regulatory chain; skeletal muscle isoform; Phosphorylase kinase alpha M subunit
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PHKA1

Image



Western blot analysis of extracts from Jurkat cells and COLO205 cells, using KPBI/2 antibody.



Immunofluorescence analysis of HeLa cells, using KP1/2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.