




YAP1 Antibody

Product Code	CSB-PA913770
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46937
Immunogen	Synthesized peptide derived from internal of Human YAP.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total YAP protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	65 kDa Yes-associated protein; YAP1; YAP65;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	YAP1

Image



Western blot analysis of extracts from HepG2 cells, COLO205 cells and HeLa cells, using YAP antibody.