




HOXA1 Antibody

Product Code	CSB-PA252418
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49639
Immunogen	Synthesized peptide derived from internal of Human HOXA1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total HOXA1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Homeobox protein Hox-A1; Hox-1F; HOXA1; HOX1F;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HOXA1

Image



Western blot analysis of extracts from COLO cells, HeLa/ cells and 293 cells, using HOXA1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.