




FLT3LG Antibody

Product Code	CSB-PA778110
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49771
Immunogen	Synthesized peptide derived from internal of Human Flt3 ligand.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total Flt3 ligand protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	SL cytokine [Precursor]; Fms-related tyrosine kinase 3 ligand; Flt3 ligand; Flt3L; FLT3LG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FLT3LG

Image



Western blot analysis of extracts from Jurkat cells, using Flt3 ligand antibody.