





FXR2 Antibody

Product Code	CSB-PA870516
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51116
Immunogen	Synthesized peptide derived from internal of Human FXR2.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total FXR2 protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Fragile X mental; retardation syndrome; related protein 2;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FXR2


Image



Western blot analysis of extracts from COLO205 cells, using FXR2 antibody.



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using FXR2 antibody.



Immunofluorescence analysis of A549 cells, using FXR2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.