




RORC Antibody

Product Code	CSB-PA439853
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51449
Immunogen	Synthesized peptide derived from internal of Human RORG.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total RORG protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Nuclear receptor ROR-gamma; Retinoid-related orphan receptor-gamma; Nuclear receptor RZR-gamma; Nuclear receptor subfamily 1 group F member 3; RORC
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RORC

Image



Western blot analysis of extracts from HuvEc cells, using RORG antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.