




ZNF75D Antibody

Product Code	CSB-PA289645
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51815
Immunogen	Synthesized peptide derived from internal of Human ZNF75.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ZNF75 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Zinc finger protein 75; Zinc finger protein 82; ZNF75; ZNF82;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZNF75D

Image



Western blot analysis of extracts from HUVEC cells, using ZNF75 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.