





EFNA5 Antibody

Product Code	CSB-PA905708
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52803
Immunogen	Synthesized peptide derived from internal of Human EFNA5.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total EFNA5 protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Ephrin-A5 [Precursor]; EPH-related receptor tyrosine kinase ligand 7; LERK-7; AL-1; EFNA5
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EFNA5


Image



Western blot analysis of extracts from HeLa cells, using EFNA5 antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using EFNA5 antibody.



Immunofluorescence analysis of A549 cells, using EFNA5 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.