




DGKQ Antibody

Product Code	CSB-PA239045
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52824
Immunogen	Synthesized peptide derived from internal of Human DGKQ.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total DGKQ protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DAG kinase theta; DAGK; DAGK4; DAGK7; DGK-theta
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DGKQ

Image



Western blot analysis of extracts from HT-29 cells, using DGKQ antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.