





CEBPG Antibody

Product Code	CSB-PA083093
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53567
Immunogen	Synthesized peptide derived from internal of Human CEBPG.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total CEBPG protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CCAAT/enhancer-binding protein gamma; C/EBP gamma; CEBPG;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CEBPG

Image



Western blot analysis of extracts from RAW264.7 cells, using CEBPG antibody.



Immunohistochemistry analysis of paraffin-embedded human heart tissue using CEBPG antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.