




EPHB4 Antibody

Product Code	CSB-PA277115
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54760
Immunogen	Synthesized peptide derived from internal of Human EPHB4.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total EPHB4 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.1.112; EC 2.7.10.1; EPB4; Ephrin type-B receptor 4 precursor; HTK
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EPHB4

Image



Western blot analysis of extracts from 293 cells and Jurkat cells, using EPHB4 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.