




HNRNPH2 Antibody

Product Code	CSB-PA061816
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55795
Immunogen	Synthesized peptide derived from internal of Human hnRNP H.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total hnRNP H protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	FTP-3; FTP3; HNRH2; HNRPH2; Heterogeneous nuclear ribonucleoprotein H
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HNRNPH2

Image



Western blot analysis of extracts from Jurkat cells, 293 cells and K562 cells, using hnRNP H antibody.