




PER3 Antibody

Product Code	CSB-PA165425
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56645
Immunogen	Synthesized peptide derived from internal of Human PER3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total PER3 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Period circadian protein homolog 3; Circadian clock protein PERIOD 3; Hper3; Cell growth-inhibiting gene 13 protein; PER3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PER3

Image



Western blot analysis of extracts from Jurkat cells, treated with insulin (0.01U/ml, 15mins), using PER3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.