





SOX8/SOX9/SOX17/SOX18 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA725661 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P57073 |
| Immunogen | Synthesized peptide derived from C-terminal of Human SOX8/9/17/18. |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse |
| Specificity | The antibody detects endogenous levels of total SOX8/9/17/18 protein. |
| Tested Applications | ELISA, WB, IHC; WB: 1:500-1:3000, IHC: 1:50-1:100 |
| Form | Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Alias | Transcription factor SOX-8/9/17/18; SOX8/9/17/18; |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | SOX8/SOX9/SOX17/SOX18 |

Image



Western blot analysis of extracts from COS7 cells, untreated or treated with serum (10%, 30mins), using SOX8/9/17/18 antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using SOX8/9/17/18 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.