





TMPRSS3 Antibody

Product Code	CSB-PA000297
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57727
Immunogen	Synthesized peptide derived from C-terminal of Human TMPRSS3.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total TMPRSS3 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Transmembrane protease; serine 3; EC 3.4.21.-; Serine protease TADG-12; Tumor-associated differentially-expressed gene 12 protein
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TMPRSS3

Image



Western blot analysis of extracts from HUVEC cells, using TMPRSS3 antibody.



Immunofluorescence analysis of HeLa cells, using TMPRSS3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.