




PKIA Antibody

Product Code	CSB-PA858680
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61925
Immunogen	Synthesized peptide derived from internal of Human IPKA.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total IPKA protein.
Tested Applications	ELISA, IHC, IF; IHC: 1:50-1:100, IF: 1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	IPKA; OKRBC1; PKI-alpha; PRKACN1; cAMP-dependent protein kinase inhibitor alpha
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PKIA

Image



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using TPD52 antibody.



Immunofluorescence analysis of HeLa cells,
using TPD52 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.